


Anti-IL6St Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P40189

Item Number: CM0020300



Introduction

Rabbit polyclonal antibody targeting human IL6ST (gp130, CD130), suitable for WB, IF/ICC, and ELISA applications. Cross-reacts with human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	IL6ST Rabbit pAb
Remarks/Alias	CD130; GP130; HIES4; IMD94; STWS2; CDW130; HIES4A; HIES4B; IL-6RB; sGP130; IL6ST
Species	Human
GeneID (Human)	3572
GeneID	3572
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 751-918 of human IL6ST (NP_002175.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40189
Protein Weight	104kDa
Shipping	Ice bag