


Anti-IL6St Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P40189

Item Number: CM0016990



Introduction

High-quality rabbit polyclonal antibody targeting human IL6ST/gp130 (CD130), validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for studying IL-6 family cytokine signaling, immune response, and inflammation.

[Learn More](#)

Parameter	Value
Product Name	IL6ST Rabbit pAb
Remarks/Alias	CD130; GP130; HIES4; IMD94; STWS2; CDW130; HIES4A; HIES4B; IL-6RB; sGP130; IL6ST
Species	Human
Gene ID (Human)	3572
Gene ID	3572
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 678-857 of human IL6ST (NP_001177910.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40189
Protein Weight	104kDa
Shipping	Ice bag