

Anti-IL7 Rabbit Monoclonal Antibody For Wb, Elisa - P13232

Item Number: CM0007463



Introduction

Recombinant rabbit monoclonal antibody targeting human IL-7. Validated for WB and ELISA. IL-7 is a key hematopoietic cytokine regulating T- and B-cell development. Protein weight: 20 kDa.

[Learn More](#)

Parameter	Value
Product Name	IL7 Rabbit mAb
Remarks/Alias	IL-7; IL7
Species	Human
GeneID	3574 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-177 of human IL7 (NP_000871.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P13232
Protein Weight	20kDa
Shipping	Ice bag