


Anti-IL8 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Fc, Elisa - P10145

Item Number: CM0006502



Introduction

Rabbit monoclonal antibody against human IL8 (Interleukin-8), validated for WB, IHC-P, IF/ICC, FC, ELISA. Targets chemokine involved in inflammation and neutrophil chemotaxis. Immunogen: recombinant protein.

[Learn More](#)

Parameter	Value
Product Name	IL8 Rabbit mAb
Remarks/Alias	IL8; NAF; GCP1; LECT; LUCT; NAP1; GCP-1; LYNAP; MDNCF; MONAP; NAP-1; SCYB8
Species	Human
Gene ID (Human)	3576
Gene ID	3576
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-99 of human IL8 (NP_000575.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, FC (intra), ELISA
Cross Reactivity	Human
SWISS	P10145
Protein Weight	11 kDa
Shipping	Ice bag