


Anti-IL8 Rabbit Polyclonal Antibody For Wb And Elisa - P10145

Item Number: CM0015279



Introduction

High-quality rabbit polyclonal antibody targeting human IL-8/CXCL8, validated for WB and ELISA. Ideal for studying neutrophil chemotaxis and inflammatory responses. Cross-reacts with human samples.

[Learn More](#)

Parameter	Specification
Product Name	IL8 Rabbit pAb
Remarks/Alias	IL8; NAF; GCP1; LECT; LUCT; NAP1; GCP-1; LYNAP; MDNCF; MONAP; NAP-1; SCYB8
Species	Human
GeneID (human)	3576
GeneID	3576
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-99 of human IL8 (NP_000575.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P10145
Protein Weight	11kDa
Shipping	Ice bag