


Anti-IL8 Rabbit Polyclonal Antibody For Wb, Elisa - P10145

Item Number: CM0020426



Introduction

Rabbit polyclonal antibody targeting human interleukin-8 (IL-8/CXCL8), validated for Western blot (WB) and ELISA. Suitable for studying chemokine-mediated inflammatory responses and neutrophil activation.

[Learn More](#)

Parameter	Value
Product Name	IL8 Rabbit pAb
Remarks/Alias	IL8; NAF; GCP1; LECT; LUCT; NAP1; GCP-1; LYNAP; MDNCF; MONAP; NAP-1; SCYB8
Species	Human
GeneID (Human)	3576
GeneID	3576
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P10145
Protein Weight	11kDa
Shipping	Ice bag