

Anti-IL9 Rabbit Monoclonal Antibody For Wb, Fc, Elisa - P15248

Item Number: CM0007697



Introduction

This recombinant rabbit monoclonal antibody targets human IL-9 (SWISS: P15248). Validated for Western blot, flow cytometry (intracellular), and ELISA, it is ideal for studying immune responses against parasites, T-cell differentiation, and mast cell function.

[Learn More](#)

Parameter	Specification
Product Name	IL9 Rabbit mAb
Remarks/Alias	P40; HP40; IL-9; IL9
Species	Human
GeneID (Human)	3578
GeneID	3578
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 19-144 of human IL9(NP_000581.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, FC (intra), ELISA
Cross Reactivity	Human
SWISS	P15248
Protein Weight	16kDa
Shipping	Ice bag