

Anti-IL9 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P15248

Item Number: CM0023614



Introduction

Rabbit polyclonal antibody against human IL9 (Interleukin-9), validated for WB, IHC-P, ELISA. Cross-reacts with rat. Suitable for immune response and cytokine research.

[Learn More](#)

Parameter	Value
Product Name	IL9 Rabbit pAb
Remarks/Alias	P40; HP40; IL-9; IL9
Species	Human
GeneID (human)	3578
GeneID	3578
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-144 of human IL9 (NP_000581.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P15248
Protein Weight	16kDa
Shipping	Ice bag