

Anti-IL9R Polyclonal Antibody (Rabbit) For Wb And Elisa - Q01113

Item Number: CM0028325



Introduction

Rabbit polyclonal antibody targeting human IL-9 receptor (IL9R, CD129). Validated for Western blotting and ELISA, with cross-reactivity to mouse and rat. Supports immunology and parasitology research.

[Learn More](#)

Parameter	Value
Product Name	IL9R Rabbit pAb
Remarks/Alias	CD129; IL-9R; IL9R
Species	Human
Gene ID (Human)	3581
Gene ID	3581
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human IL9R (XP_011529456.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01113
Protein Weight	57kDa
Shipping	Ice bag