

Anti-Ilf2 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q12905

Item Number: CM0010930

Introduction

Rabbit monoclonal antibody to human ILF2. Validated for WB, IHC-P, ELISA with broad reactivity (human, mouse, rat). Ideal for transcription and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	ILF2 Rabbit mAb
Remarks/Alias	NF45; PRO3063; ILF2
Species	Human
GeneID	3608
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human ILF2 (Q12905)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12905
Protein Weight	43kDa
Shipping	Ice bag