

## Anti-ILf3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q12906

Item Number: CM0004174



### Introduction

Recombinant rabbit monoclonal antibody against human ILF3 (UniProt Q12906). Validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Detects a 95kDa band. Ideal for investigating circRNA biogenesis, transcriptional regulation, and innate antiviral responses.

[Learn More](#)

| Specification    | Value                                                                                                                                                                             |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | ILF3 Rabbit mAb                                                                                                                                                                   |
| Remarks/Alias    | CBTF; DRBF; MMP4; MPP4; NF90; NFAR; NF110; NF90a; NF90b; NF90c; NFAR2; TCP80; DRBP76; NF110b; NFAR-1; NFAR-2; NFAR90; TCP110; MPP4110; NF90ctv; NFAR110; MPHOSPH4; NF-AT-90; ILF3 |
| Species          | Human                                                                                                                                                                             |
| GeneID (Human)   | 3609                                                                                                                                                                              |
| GeneID           | 3609                                                                                                                                                                              |
| Immunogen        | Synthetic peptide corresponding to amino acids 1-100 of human ILF3 (Q12906)                                                                                                       |
| Source           | Rabbit                                                                                                                                                                            |
| Category         | Monoclonal Antibodies                                                                                                                                                             |
| Application      | WB, IHC-P, ELISA                                                                                                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                                                                 |
| SWISS            | Q12906                                                                                                                                                                            |
| Protein Weight   | 95kDa                                                                                                                                                                             |
| Shipping         | Ice bag                                                                                                                                                                           |