

Anti-Ilk Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13418

Item Number: CM0012869



Introduction

ILK Rabbit pAb is a polyclonal antibody targeting scaffold protein ILK, validated for WB, IHC-P, and ELISA. It cross-reacts with human, mouse, and rat ILK, supporting studies in cell adhesion, migration, and focal adhesion dynamics.

[Learn More](#)

Parameter	Value
Product Name	ILK Rabbit pAb
Remarks/Alias	P59; ILK-1; ILK-2; p59ILK; HEL-S-28; ILK
Species	Human
GeneID (Human)	3611
GeneID	3611
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human ILK (NP_001014795.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13418
Protein Weight	51kDa
Shipping	Ice bag