

Anti-Ilkap Polyclonal Antibody For Wb And Elisa - Q9H0C8

Item Number: CM0013452

Introduction

Rabbit polyclonal antibody against ILKAP (PP2C-DELTA), validated for Western blot and ELISA, cross-reacts with human and mouse. Ideal for studying cell adhesion, growth factor signaling, and cell cycle phosphatase activity.

[Learn More](#)

Parameter	Value
Product Name	ILKAP Rabbit pAb
Remarks/Alias	PP2CD; PPM1O; ILKAP2; ILKAP3; PP2C-DELTA; ILKAP
Species	Human
GeneID (Human)	80895
GeneID	80895
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-392 of human ILKAP (NP_110395.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9H0C8
Protein Weight	43kDa
Shipping	Ice bag