


Anti-IlT5/Cd85A (Lilrb3) Monoclonal Antibody () For Flow Cytometry - O75022

Item Number: CM0011595



Introduction

-conjugated rabbit monoclonal antibody against human ILT5/CD85a (LILRB3) for flow cytometry. Validated in human samples; detects the 69-71 kDa inhibitory receptor on monocytes and B-cells.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ILT5/CD85a mAb
Remarks/Alias	LILRB3; CD85A; HL9; ILT-5; ILT5; LILRA6; LIR-3; LIR3; PIR-B; PIRB; leukocyte immunoglobulin like receptor B3
Species	Human
Gene ID (Human)	11025
Gene ID	11025
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-418 of human ILT5/CD85a (NP_006855.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O75022
Protein Weight	69kDa/71kDa
Shipping	Ice bag