

Anti-Immp2L Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q96T52

Item Number: CM0019354

Introduction

Rabbit polyclonal antibody against human IMMP2L (IMP2-like protein), a mitochondrial inner membrane protease involved in protein import and apoptosis. Validated for WB, IF/ICC, ELISA; reacts with human and mouse samples.

[Learn More](#)

Parameter	Value
Product Name	IMMP2L Rabbit pAb
Remarks/Alias	IMMP2L; IMMP2L-IT1; IMP2; IMP2-LIKE; mitochondrial inner membrane protease subunit 2
Species	Human
GeneID (Human)	83943
GeneID	83943
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 38-175 of human IMMP2L(NP_001231535.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q96T52
Protein Weight	12kDa/19kDa
Shipping	Ice bag