

## Anti-Immt Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc - Q16891

Item Number: CM0027135

### Introduction

Rabbit polyclonal antibody targeting human IMMT (Mitofilin, MIC60) for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with human, mouse, rat. Ideal for mitochondrial inner membrane architecture and cristae morphology research.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | IMMT Rabbit pAb                                                                                                    |
| Remarks/Alias    | HMP; P87; P89; PIG4; Mic60; PIG52; MINOS2; P87/89; MICOS60; IMMT                                                   |
| Species          | Human                                                                                                              |
| Gene ID          | 10989                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 120-270 of human IMMT (NP_006830.2). |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                           |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | Q16891                                                                                                             |
| Protein Weight   | 84kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |