


Anti-Impa1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P29218

Item Number: CM0026471



Introduction

IMPA1 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA in human, mouse, rat. Targets inositol monophosphatase 1 (IMPA1), a key enzyme in phosphatidylinositol signaling and lithium-sensitive pathway. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	IMPA1 Rabbit pAb
Remarks/Alias	IMP; IMPA; MRT59; IMPA1
Species	Human
Gene ID	3612
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-277 of human IMPA1 (NP_005527.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P29218
Protein Weight	30kDa
Shipping	Ice bag