


## Anti-Impdh2 Polyclonal Antibody For Wb, If/Icc, Elisa - P12268

Item Number: CM0025319



### Introduction

IMPDH2 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat IMPDH2 (~56 kDa), a key enzyme in guanine nucleotide synthesis. Ideal for cancer research and nucleotide metabolism studies.

[Learn More](#)

| Parameter        | Value                                                                                                               |
|------------------|---------------------------------------------------------------------------------------------------------------------|
| Product Name     | IMPDH2 Rabbit pAb                                                                                                   |
| Remarks/Alias    | IMPD2; IMPDH-II; IMPDH2                                                                                             |
| Species          | Human                                                                                                               |
| GeneID           | 3615                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 215-514 of human IMPDH2 (NP_000875.2) |
| Source           | Rabbit                                                                                                              |
| Category         | Polyclonal Antibodies                                                                                               |
| Application      | WB, IF/ICC, ELISA                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                   |
| SWISS            | P12268                                                                                                              |
| Protein Weight   | 56kDa                                                                                                               |
| Shipping         | Ice bag                                                                                                             |