


Anti-Impdh2 Polyclonal Antibody For Wb, Ihc-P, Elisa - P12268

Item Number: CM0017531



Introduction

IMPDH2 rabbit polyclonal antibody validated for Western blot, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Target inosine-5'-monophosphate dehydrogenase 2, a key enzyme in de novo guanine nucleotide synthesis (SWISS P12268).

[Learn More](#)

Product Name	IMPDH2 Rabbit pAb
Remarks/Alias	IMPD2; IMPDH-II; IMPDH2
Species	Human
GeneID (Human)	3615
GeneID	3615
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 215-514 of human IMPDH2 (NP_000875.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P12268
Protein Weight	56kDa
Shipping	Ice bag