


Anti-Importin 9 / Ranbp9 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q96P70

Item Number: CM0003619



Introduction

Rabbit monoclonal antibody targeting Importin 9 (RANBP9), validated for Western blot, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for studying nuclear transport and histone chaperone functions.

[Learn More](#)

Parameter	Value
Product Name	Importin 9 / RANBP9 Rabbit mAb
Remarks/Alias	Imp9; Importin 9 / RANBP9
Species	Human
Gene ID (Human)	55705
Gene ID	55705
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human Importin 9 / RANBP9 (Q96P70)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96P70
Protein Weight	116 kDa
Shipping	Ice bag