

Anti-Importin Alpha 5/Kpna1 Rabbit Mab - P52294

Item Number: CM0007261



Introduction

Monoclonal antibody against human Importin Alpha 5 (KPNA1) validated for WB, IF/ICC, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for studying nucleocytoplasmic transport mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Importin Alpha 5/KPNA1 Rabbit mAb
Remarks/Alias	RCH2; SRP1; IPOA5; NPI-1; Importin Alpha 5/KPNA1
Species	Human
GeneID (Human)	3836
GeneID	3836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human Importin Alpha 5/KPNA1 (NP_002255.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52294
Protein Weight	60kDa
Shipping	Ice bag