


Anti-Inadl Rabbit Polyclonal Antibody For Wb, Elisa - Q8Ni35

Item Number: CM0024778



Introduction

INADL Rabbit pAb is a rabbit polyclonal antibody raised against human INADL using a synthetic peptide immunogen. It is recommended for WB and ELISA, recognizing human and mouse INADL, a tight junction scaffolding protein.

[Learn More](#)

Parameter	Value
Product Name	INADL Rabbit pAb
Remarks/Alias	Cipp; INADL; hINADL; InaD-like
Species	Human
Gene ID (Human)	10207
Gene ID	10207
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q8Ni35
Protein Weight	196kDa
Shipping	Ice bag