


Anti-Incenp Polyclonal Antibody For Wb And Elisa - Q9Nqs7

Item Number: CM0018338



Introduction

Rabbit polyclonal antibody against human INCENP, validated for Western blot and ELISA. Suitable for detecting INCENP in human and rat samples. Ideal for mitosis and chromosome segregation research.

[Learn More](#)

Parameter	Specification
Product Name	INCENP Rabbit pAb
Remarks/Alias	INCENP
Species	Human
GeneID (human)	3619
GeneID	3619
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 785-914 of human INCENP (NP_064623.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9NQS7
Protein Weight	105kDa
Shipping	Ice bag