


Anti-Indolethylamine N-Methyltransferase (Temt) Rabbit Monoclonal Antibody For Wb, Elisa - O95050

Item Number: CM0002510



Introduction

High-quality rabbit monoclonal antibody targeting human Indolethylamine N-methyltransferase (TEM), suitable for Western blot and ELISA. Detects a band at ~29 kDa; cross-reacts with mouse. Ideal for methylation and detoxification research.

[Learn More](#)

Parameter	Value
Product Name	TEM Rabbit mAb
Remarks/Alias	TEM
Species	Human
Gene ID (Human)	11185
Gene ID	11185
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human TEM (O95050)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O95050
Protein Weight	29kDa
Shipping	Ice bag