


Anti-Ing1 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Uk53

Item Number: CM0027577



Introduction

Rabbit polyclonal antibody targeting human ING1, validated for WB and ELISA. ING1 is a tumor suppressor cooperating with p53 in cell growth regulation. Predicted MW: 47kDa.

[Learn More](#)

Parameter	Value
Product Name	ING1 Rabbit pAb
Remarks/Alias	p33; p47; p33ING1; p24ING1c; p33ING1b; p47ING1a; ING1
Species	Human
GeneID (Human)	3621
GeneID	3621
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 80-180 of human ING1 (NP_937862.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Uk53
Protein Weight	47kDa
Shipping	Ice bag