


## Anti-Ing4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Unl4

Item Number: CM0030974



### Introduction

Rabbit polyclonal antibody targeting human ING4 (p29ING4), validated for WB, IHC-P, and ELISA. Reacts with human and mouse samples; ideal for tumor suppression, chromatin, and apoptosis research.

[Learn More](#)

| Parameter        | Value                 |
|------------------|-----------------------|
| Product Name     | ING4 Rabbit pAb       |
| Remarks/Alias    | my036; p29ING4; ING4  |
| Species          | Human                 |
| GeneID (Human)   | 51147                 |
| GeneID           | 51147                 |
| Immunogen        | `Recombinant protein  |
| Source           | Rabbit                |
| Category         | Polyclonal Antibodies |
| Application      | WB, IHC-P, ELISA      |
| Cross Reactivity | Human, Mouse          |
| SWISS            | Q9UNL4                |
| Protein Weight   | 29kDa                 |
| Shipping         | Ice bag               |