


Anti-Inhbb Rabbit Polyclonal Antibody For Wb, Elisa - P09529

Item Number: CM0030145



Introduction

Rabbit polyclonal antibody targeting human INHBB (Inhibin beta B chain). Validated for Western blot and ELISA. Ideal for reproductive biology, stem cell, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	INHBB Rabbit pAb
Remarks/Alias	INHBB
Species	Human
GeneID (human)	3625
GeneID	3625
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 293-407 of human INHBB (NP_002184.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P09529
Protein Weight	45kDa
Shipping	Ice bag