


Anti-Inhibin Beta A (Inhba) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P08476

Item Number: CM0030513



Introduction

Rabbit polyclonal antibody against human Inhibin beta A (INHBA). Validated for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for TGF-beta signaling, reproductive biology, and developmental research.

[Learn More](#)

Parameter	Value
Product Name	Inhibin beta A (INHBA) Rabbit pAb
Remarks/Alias	EDF; FRP; Inhibin beta A (INHBA)
Species	Human
GeneID	3624
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-426 of human Inhibin beta A (INHBA) (NP_002183.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08476
Protein Weight	47kDa
Shipping	Ice bag