


Anti-Inhibin Beta A Rabbit Polyclonal Antibody For Wb, Elisa - P08476

Item Number: CM0029888



Introduction

Rabbit polyclonal antibody targeting human Inhibin beta A (INHBA), cross-reactive with mouse. Suitable for Western blot and ELISA, supporting research into reproductive biology, stem cell maintenance, and fibrosis. Provided by CamelBio.

[Learn More](#)

Parameter	Value
Product Name	Inhibin beta A Rabbit pAb
Remarks/Alias	EDF; FRP; Inhibin beta A (INHBA)
Species	Human
GeneID (human)	3624
GeneID	3624
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Inhibin beta A (NP_002183.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08476
Protein Weight	47kDa
Shipping	Ice bag