

Anti-Inos Monoclonal Antibody For Wb, If/Icc, Elisa - P35228

Item Number: CM0010644



Introduction

Rabbit monoclonal antibody targeting human iNOS (NOS2), suitable for Western blot, immunofluorescence, and ELISA. Cross-reacts with mouse iNOS. Ideal for studying inflammation, immune response, and nitric oxide signaling.

[Learn More](#)

Parameter	Value
Product Name	iNOS Rabbit mAb
Remarks/Alias	NOS; INOS; NOS2A; HEP-NOS
Species	Human
GeneID (Human)	4843
GeneID	4843
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P35228
Protein Weight	131kDa
Shipping	Ice bag