

Anti-Inos Rabbit Monoclonal Antibody For Wb - P35228

Item Number: CM0012183



Introduction

Rabbit monoclonal antibody against human iNOS (P35228), suitable for Western blot, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Produced using recombinant protein immunogen (aa 2-200).

[Learn More](#)

Parameter	Value
Product Name	iNOS Rabbit mAb
Remarks/Alias	NOS; INOS; NOS2A; HEP-NOS; iNOS
Species	Human
GeneID (Human)	4843
GeneID	4843
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2-200 of human iNOS (P35228)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P35228
Protein Weight	131 kDa
Shipping	Ice bag