

Anti-Inos Rabbit Polyclonal Antibody For Wb And Elisa - P35228

Item Number: CM0016394



Introduction

Rabbit polyclonal antibody against human iNOS (inducible nitric oxide synthase) for WB and ELISA; cross-reacts with human and mouse. Suitable for inflammation, macrophage biology, and NO signaling research.

[Learn More](#)

Parameter	Value
Product Name	iNOS Rabbit pAb
Remarks/Alias	NOS; INOS; NOS2A; HEP-NOS; iNOS
Species	Human
GeneID (Human)	4843
GeneID	4843
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 974-1143 of human iNOS (NP_000616.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P35228
Protein Weight	131kDa
Shipping	Ice bag