

Anti-Inpp1 Rabbit Polyclonal Antibody For Wb And Elisa - P49441

Item Number: CM0013128

Introduction

Rabbit polyclonal antibody targeting human INPP1 (inositol polyphosphate 1-phosphatase), suitable for Western blot and ELISA. Validated for human reactivity; immunogen covers residues 1-180.

[Learn More](#)

Parameter	Value
Product Name	INPP1 Rabbit pAb
Remarks/Alias	INPP1
Species	Human
GeneID (Human)	3628
GeneID	3628
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human INPP1 (NP_002185.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P49441
Protein Weight	44kDa
Shipping	Ice bag