

Anti-Inpp4A Polyclonal Antibody For Wb, Elisa - Q96Pe3

Item Number: CM0029878

Introduction

Rabbit polyclonal antibody targeting human INPP4A (Q96PE3), suitable for WB and ELISA. Cross-reacts with mouse and rat. Detects the N-terminal region to study phosphoinositide signaling and neuroprotection.

[Learn More](#)

Parameter	Value
Product Name	INPP4A Rabbit pAb
Remarks/Alias	INPP4; TVAS1; INPP4A
Species	Human
GeneID (Human)	3631
GeneID	3631
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human INPP4A (NP_004018.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96PE3
Protein Weight	110kDa
Shipping	Ice bag