

Anti-Inpp4B Rabbit Polyclonal Antibody For Wb And Elisa - O15327

Item Number: CM0022071

Introduction

Rabbit polyclonal antibody against human INPP4B, validated for WB and ELISA. Cross-reacts with mouse and rat. Suitable for studying lipid phosphatase function and PI3K-AKT signaling.

[Learn More](#)

Parameter	Value
Product Name	INPP4B Rabbit pAb
Remarks/Alias	INPP4B
Species	Human
Gene ID (human)	8821
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 380-600 of human INPP4B (NP_003857.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15327
Protein Weight	105kDa
Shipping	Ice bag