

## Anti-Inpp5A Polyclonal Antibody For Western Blot And Elisa - Q14642

Item Number: CM0025293

### Introduction

High-quality rabbit polyclonal antibody targeting human INPP5A (5PTase) for inositol phosphate signaling studies. Validated for WB and ELISA with mouse and rat cross-reactivity. Ideal for heart, brain, and muscle tissue analysis.

[Learn More](#)

| Parameter        | Value                                                                                                                                   |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | INPP5A Rabbit pAb                                                                                                                       |
| Remarks/Alias    | 5PTASE; INPP5A                                                                                                                          |
| Species          | Human                                                                                                                                   |
| Gene ID (Human)  | 3632                                                                                                                                    |
| Gene ID          | 3632                                                                                                                                    |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 10-200 of human INPP5A (NP_005530.3). |
| Source           | Rabbit                                                                                                                                  |
| Category         | Polyclonal Antibodies                                                                                                                   |
| Application      | WB, ELISA                                                                                                                               |
| Cross Reactivity | Mouse, Rat                                                                                                                              |
| SWISS            | Q14642                                                                                                                                  |
| Protein Weight   | 48kDa                                                                                                                                   |
| Shipping         | Ice bag                                                                                                                                 |