

Anti-Inpp5E Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Nrr6

Item Number: CM0023841



Introduction

High-quality anti-INPP5E rabbit polyclonal antibody validated for WB, IF/ICC, ELISA, cross-reactive with human, mouse, rat. Ideal for studying phosphoinositide signaling and ciliopathy research.

[Learn More](#)

Parameter	Value
Product Name	INPP5E Rabbit pAb
Remarks/Alias	CPD4; CORS1; JBTS1; MORMS; PPI5PIV; pharbin; INPP5E
Species	Human
GeneID (Human)	56623
GeneID	56623
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 510-630 of human INPP5E (NP_063945.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRR6
Protein Weight	70kDa
Shipping	Ice bag