

Anti-Inpp5J Polyclonal Antibody For Wb, If/Icc, Elisa - Q15735

Item Number: CM0027831

Introduction

Rabbit polyclonal antibody against INPP5J (PIPP, INPP5) for WB, IF/ICC, ELISA. Validated in human, mouse, rat. Target weight 107kDa. SWISS Q15735. Ideal for phosphoinositide signaling research.

[Learn More](#)

Parameter	Value
Product Name	INPP5J Rabbit pAb
Remarks/Alias	PIPP, INPP5, PIB5PA, INPP5J
Species	Human
GeneID (Human)	27124
GeneID	27124
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 60-360 of human INPP5J (NP_001002837.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15735
Protein Weight	107kDa
Shipping	Ice bag