

Anti-Inpp5K Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bt40

Item Number: CM0029693

Introduction

Anti-INPP5K Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA in Human, Mouse, Rat. Ideal for studying inositol phosphate signaling, actin cytoskeleton regulation, and insulin-mediated glucose uptake.

[Learn More](#)

Parameter	Value
Product Name	INPP5K Rabbit pAb
Remarks/Alias	PPS; SKIP; MDCC4ID; INPP5K
Species	Human
GeneID	51763
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 149-448 of human INPP5K (NP_057616.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BT40
Protein Weight	51kDa
Shipping	Ice bag