

Anti-Insig2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Y5U4

Item Number: CM0021114



Introduction

Rabbit polyclonal antibody targeting human INSIG2, validated for WB, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat. INSIG2 is a key oxysterol-binding protein regulating cholesterol synthesis via SREBP retention and HMGCR degradation.

[Learn More](#)

Parameter	Value
Product Name	INSIG2 Rabbit pAb
Remarks/Alias	INSIG-2; INSIG2
Species	Human
GeneID (Human)	51141
GeneID	51141
Immunogen	Synthetic peptide — A synthetic peptide corresponding to a sequence within amino acids 50-150 of human INSIG2 (NP_057217.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5U4
Protein Weight	25kDa
Shipping	Ice bag