

Anti-Insr/Igf1R Monoclonal Antibody For Wb And Elisa - P06213 / P08069

Item Number: CM0007720



Introduction

Anti-INSR/IGF1R Rabbit mAb for WB and ELISA, targeting human, mouse, rat insulin/IGF-1 receptors (156 kDa). Ideal for signal transduction, metabolism, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	INSR/IGF1R Rabbit mAb
Remarks/Alias	HHF5; CD220; IGFR; CD221; IGFR; JTK13; INSR/IGF1R
Species	Human
GeneID (Human)	3643/3480
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1100-1200 of human INSR/IGF1R. (NP_000199.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06213/P08069
Protein Weight	156kDa
Shipping	Ice bag