

Anti-Insulin Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P01308

Item Number: CM0002847



Introduction

Rabbit monoclonal antibody against human Insulin, validated for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for diabetes and metabolic disorder studies.

[Learn More](#)

Parameter	Value
Product Name	Insulin Rabbit mAb
Remarks/Alias	IDDM; ILPR; IRDN; IDDM1; IDDM2; PNDM4; MODY10; Insulin
Species	Human
Gene ID (Human)	3630
Gene ID	3630
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 25-110 of human Insulin (P01308).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01308
Protein Weight	12kDa
Shipping	Ice bag