

Anti-Insulin Receptor Polyclonal Antibody For Wb, Elisa - P06213

Item Number: CM0020707



Introduction

Rabbit polyclonal antibody against human Insulin Receptor (P06213) validated for Western Blot and ELISA. Cross-reacts with mouse and rat. Suitable for insulin signaling and metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	Insulin Receptor Rabbit pAb
Remarks/Alias	HHF5; CD220; Insulin Receptor
Species	Human
GeneID (Human)	3643
GeneID	3643
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 475-570 of human Insulin Receptor (NP_000199.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human,Mouse,Rat
SWISS	P06213
Protein Weight	156kDa
Shipping	Ice bag