

Anti-Insulin Receptor Rabbit Monoclonal Antibody For Wb, Elisa - P06213

Item Number: CM0002473



Introduction

Recombinant rabbit monoclonal antibody against human insulin receptor (P06213), validated for WB and ELISA, cross-reacts with human, mouse, and rat. Ideal for insulin signaling and metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	Insulin Receptor Rabbit mAb
Remarks/Alias	HHF5; CD220; Insulin Receptor
Species	Human
Gene ID (Human)	3643
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 900-1000 of human Insulin Receptor (P06213).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06213
Protein Weight	156kDa
Shipping	Ice bag