

Anti-Insulin Receptor Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P06213

Item Number: CM0020304



Introduction

Rabbit polyclonal antibody against human insulin receptor (INSR), validated for Western blot (WB), immunoprecipitation (IP), and ELISA. Cross-reacts with human, mouse, and rat. Ideal for studying insulin signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Insulin Receptor Rabbit pAb
Remarks/Alias	HHF5; CD220; Insulin Receptor
Species	Human
GeneID (Human)	3643
GeneID	3643
Immunogen	Recombinant protein or recombinant fusion protein containing a sequence corresponding to amino acids 1281-1382 of human INSR (NP_000199.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06213
Protein Weight	156kDa
Shipping	Ice bag