


Anti-Integrin Alpha 1 (Itga1/Cd49A) Rabbit Monoclonal Antibody For Wb And Elisa - P56199

Item Number: CM0006361



Introduction

Rabbit monoclonal antibody targeting human Integrin alpha-1 (ITGA1/CD49a), validated for WB and ELISA. Laminin/collagen receptor involved in cell adhesion. Ideal for research and diagnostic applications.

[Learn More](#)

Parameter	Value
Product Name	Integrin alpha 1 (ITGA1/CD49a) Rabbit mAb
Remarks/Alias	VLA1; CD49a; Integrin alpha 1 (ITGA1/CD49a)
Species	Human
Gene ID (Human)	3672
Gene ID	3672
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-1141 of human Integrin alpha 1 (ITGA1/CD49a) (NP_852478.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P56199
Protein Weight	131 kDa
Shipping	Ice bag