

Anti-Integrin Alpha 1/Cd49A Monoclonal Antibody For Wb And If - Q3V3R4

Item Number: CM0010615

Introduction

Rabbit monoclonal antibody against mouse Integrin alpha 1 (CD49a/VLA-1) validated for WB, IF/ICC, IF-P, ELISA. Recognizes mouse and rat. Research use in cell adhesion and integrin signaling.

[Learn More](#)

Parameter	Value
Product Name	Integrin alpha 1/CD49a Rabbit mAb
Remarks/Alias	Vla1; CD49A; E130012M19Rik
Species	Mouse
GeneID (human)	3672
GeneID	109700
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-1141 of mouse Integrin alpha 1/CD49a (NP_001028400.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q3V3R4
Protein Weight	131 kDa
Shipping	Ice bag