


Anti-Integrin Alpha 1/Cd49A Rabbit Polymab® For Wb, Elisa - P56199 / Q3V3R4

Item Number: CM0008941



Introduction

High-quality anti-Integrin alpha 1/CD49a Rabbit PolymAb® monoclonal antibody validated for WB and ELISA, reactive with human, mouse, and rat. Ideal for studying laminin/collagen receptor signaling in cell adhesion and growth regulation.

[Learn More](#)

Parameter	Value
Product Name	Integrin alpha 1/CD49a Rabbit PolymAb®
Remarks/Alias	VLA1; CD49a
Species	Human/Mouse
Gene ID (Human)	3672
Gene ID (Mouse)	109700
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-1141 of human/mouse Integrin alpha 1/CD49a (NP_852478.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56199/Q3V3R4
Protein Weight	131kDa
Shipping	Ice bag