

Anti-Integrin Alpha 3/Cd49C Monoclonal Antibody For IHC-P, Fc, Elisa - P26006

Item Number: CM0011416

Introduction

Rabbit monoclonal antibody against human Integrin alpha 3/CD49c, suitable for IHC-P, FC, ELISA. Recognizes a 117kDa protein involved in cell adhesion, invasion, and microbial infection. Ideal for research use.

[Learn More](#)

Parameter	Value
Product Name	Integrin alpha 3/CD49c Rabbit mAb
Remarks/Alias	JEB7; VL3A; CD49C; FRP-2; GAPB3; ILNEB; MSK18; VCA-2; VLA3a; GAP-B3
Species	Human
GeneID (Human)	3675
GeneID	3675
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 466-620 of human Integrin alpha 3/CD49c (NP_002195.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, FC, ELISA
Cross Reactivity	Human
SWISS	P26006
Protein Weight	117kDa
Shipping	Ice bag