


Anti-Integrin Alpha 4 (Itga4/Cd49D) Monoclonal Antibody For Wb, Ihc-P, Elisa - P13612

Item Number: CM0012740



Introduction

Rabbit monoclonal antibody targeting human ITGA4 (CD49d), validated for WB, IHC-P, ELISA. Cross-reacts with mouse. A key reagent for studying integrin-mediated adhesion and leukocyte trafficking pathways.

[Learn More](#)

Parameter	Value
Product Name	Integrin alpha 4 (ITGA4/CD49d) Rabbit mAb
Remarks/Alias	IA4; CD49D; Integrin alpha 4 (ITGA4/CD49d)
Species	Human
Gene ID (Human)	3676
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 933-1032 of human Integrin alpha 4 (ITGA4/CD49d) (P13612)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P13612
Protein Weight	115kDa
Shipping	Ice bag